RQflex[®] plus 10

Operating manual

en

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de

Mode d'emploi

fr

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es

Contents


1	Description of the instrument	4
	Display and operation keys	4
	Technical data	6
	Accessories	6
2	Setting up the instrument	7
	Insertion of batteries	7
	Setting time and date	8
	Initial calibration	9
3	Calibration	10
4	Measurement	12
	Measurement procedure A	12
	Measurement procedure B	14
	Measurement procedure C	16
5	Handling measurement results	20
	Storing measurement results	20
	Displaying stored data	20
	Data transfer to a PC	22
	Deleting data	22
6	Maintenance of the instrument	24
	Handling	24
	Cleaning the strip adapter	25
	Cleaning the cell adapter	26
	Recalibration	26
7	Trouble-shooting	28
	Error messages, problems, solutions	28
8	Service/Warranty	30
	Service department, warranty conditions	30

1 Description of the instrument

Your RQflex® plus 10 is a versatile, precise instrument. It is part of the Reflectoquant® system with the components:

- instrument RQflex® plus 10
- Reflectoquant® tests
- test- and batch-specific bar-code strip

According to the principle of reflectometry (remission photometry), diffuse reflected light from the test strip or after transmission through a cell is measured. Just as in classical photometry, the difference in intensity of emitted and reflected light allows a quantitative determination of the concentration of specific analytes.

When the instrument is switched on, all elements of the display are shown for about one second. In this time the instrument runs a self-test. Ensure that all functions, especially the multifunctional digit display, are shown correctly.

en Display and operation keys (buttons)

Operation keys (buttons) 1

Display for messages and information 2

Interface 3

For RQdata only (see page 22)

Do not use to connect an external power supply!

Lid 4


Test-method memory 5

Up to five test methods, e.g. nitrate, pH value, can be stored simultaneously.

Test methods are loaded by means of the bar-code strip.

Bar-code scanner 6

All information necessary to run the instrument is communicated via the bar-code strip. The bar code contains information for wavelength correction and for the batch-specific calibration curve.

Time 7

Change of battery (LOW BAT) 8


Date 9

Symbols for measurement procedures 10

Test-method memory 11

Multifunctional display 12

Instrument, tests, and bar-code strips are tuned to each other. Ensure that the reference numbers (code) on the display, the bar code, the tube containing the test strips, and the test package matches.

Technical data

Dimension:	19 x 8 x 2 cm
Weight:	275 g
Memory:	5 test methods, 50 measurement results
Interface:	yes
Light source:	4 LEDs, green/red, double optics
Power source:	4 x 1.5-V batteries
Display:	LCD
System diagnosis:	yes
Measurement range:	4 – 90 % rel. remission
Reflection area:	4 x 6 mm
Resolution:	0.1 % rel. remission
Photom. accuracy:	0.5 % rel. remission (measurement to measurement, instrument to instrument)
Operating temperature:	5 – 40 °C for ideal measurements
Operating humidity:	below 90 % für for ideal measurements

Accessories

Reflectoquant® Software RQdata	1.16998.0001	for Windows incl. interface (RS 232, USB) for transferring data from RQflex® reflectometers to PC
Reflectoquant® Strip adapter	1.16953.0001	Replacement part for RQflex® and RQflex® plus
Reflectoquant® Recalibration Set	1.16954.0001	Replacement part for RQflex®
Reflectoquant® Check Set RQcheck	1.16957.0001	Accessory for RQflex® and RQflex® plus
Reflectoquant® Cell adapter	1.16729.0001	Accessory for RQflex® plus
Reflectoquant® Empty cells	1.16727.0001	Accessory for RQflex® plus

2 Setting up the instrument

Insert the four supplied batteries into the instrument.

Insertion of batteries

When?


Before first use or when the LOW BAT sign is shown in the display.

How?

1. Remove the lid of the battery compartment on the back of the instrument by pressing carefully in the indicated direction.
2. Insert the batteries into the compartment, heeding the + and – pole signs.
3. Close the battery compartment.

Important

When the LOW BAT sign appears the first time, at least a further 20 measurements can still be performed. Afterwards the instrument cannot be switched on. In order to ensure that stored measurement results are saved, the battery should be exchanged as soon as possible after LOW BAT first appears in the display. Measurement results can be irretrievably lost if the current is interrupted for more than 2 minutes, i.e. the batteries must be replaced within the space of 2 minutes.


Setting time and date

When?

1. With a new instrument.
2. With the change of summer/winter time.
3. Check after replacement of batteries.



How?

1. Switch off the instrument by pressing the ON/OFF button.
2. Press the MEM button and simultaneously press the ON/OFF key. The time and date display blinks.
3. To select the international or American date press the START button.
4. Press the TEST button to access hour, minute, day, and month.
5. Press the START button to adjust the time (hour, minute) and the date (day, month).
6. End the procedure by pressing the MEM or the ON/OFF key. The set values are automatically saved, the clock activated.

Initial calibration

1. Keep the instrument at room temperature for at least 30 minutes.
2. Switch the instrument on and insert the bar code of the recalibration set (supplied) into the bar-code scanner. CAL appears in the display.
3. Insert the calibration strip of the recalibration set into the strip adapter as shown in the diagram, and ensure that it is correctly positioned. Take utmost care to ensure that the calibration strip is clean.
4. Press the START key. The instrument is recalibrated, the CAL message disappears from the display, and the instrument switches off automatically.
5. Store the recalibration set in the supplied packaging materials, taking care to protect it from light.

Important

When the batteries are low (display ) the instrument will not accept the bar-code strip for calibration. Change the batteries first. If the calibration strip is not inserted into the measurement chamber, an acoustic signal sounds after the START button is pressed. The blinking CAL symbol in the display prompts you to insert the calibration strip.


3 Calibration

When?

1. With a new instrument.
2. With a new batch of Reflectoquant® test.
The supplied bar code is valid only for the batch with which it is enclosed.
3. When scanning in a new test method of a new test.

How?

1. Remove the specific bar-code strip from the Reflectoquant® pack. Certain tests require two bar-code strips.
2. Switch on the instrument by pressing the ON/OFF button and open the lid. The display of F:50 indicates that 50 memory slots are available.
3. Press the TEST button until the arrow points to the corresponding test-method memory where you wish to store the calibration.
4. Insert the bar code all the way into the bar-code scanner in the direction of the arrow in a smooth motion from left to right, and then remove it again.
5. The calibration process is successfully finished when the left-hand three digits of the batch number on the bar-code strip are displayed. In addition an acoustic signal sounds. The small arrow in the display indicates the test-method memory place (1 – 5) where the test data are stored.
6. Repeat this procedure until the acoustic signal sounds. Entry is now complete.
7. Return the bar-code strip to the Reflectoquant® pack. (Store carefully, but not in the test-strip tube.)

Important

Overwriting a test method

- An already calibrated test method stored in the memory is overwritten when a new bar code is scanned in.
- All measurement data stored for a given method is lost when it is overwritten.
- The blinking MEM in the display indicates that measurement data are stored. If these data are still required, they can be individually called up by pressing the MEM key and correspondingly noted.

Test methods that requires two bar-code strips

- After scanning in the first bar-code strip, test methods needing two bar-code strips show the three-digit reference number and the symbol --- in the display alternately.
- As soon as the second bar-code strip has been scanned in, only the batch reference number is shown.
- The sequence of scanning is of no importance.

4 Measurement

There are two distinct measurement procedures for test strips and one for cell tests. The both procedures for the test strips are very similar and differ only in the first steps. Each measurement procedure is test-specific. The procedure is transferred to the instrument via the bar-code strip. All steps involved to obtain the result are displayed on the screen. Please also refer to the respective Reflectoquant® package insert!

Measurement procedure A


In this procedure, only the reaction time of the Reflectoquant® test has to be considered.

How?

1. Switch on the instrument by pressing the ON/OFF button.
2. Press the TEST button until the arrow indicates the desired test-method memory.
3. Compare the first three digits of the batch number of the Reflectoquant® pack. This reference number must be the same as the one displayed (otherwise recalibrate, see page 10).
4. Press the START button. The reaction time is shown on the display.
5. Immerse the test strip into your sample as described in the Reflectoquant® package insert (Reflectoquant® pack) and at the same time press the START button (A). Carefully allow excess liquid to run off via the long edge of the strip onto an absorbent paper towel. The countdown function of the instrument is activated, the remaining reaction time is shown (countdown) (B).

6. An acoustic signal and the blinking strip five seconds before the end of the reaction time requests you to insert the strip all the way into the strip adaptor. The strip can also be inserted before the five seconds have expired (see the corresponding package insert).
7. The measurement result is shown in the display in the units for the respective test, e.g. mg/l, and is automatically stored.

2

5B


4


6

5A

7







Measurement procedure B

Besides the reaction time of the Reflectoquant® tests this procedure also considers other waiting times.

How?

1. Switch on the instrument by pressing the ON/OFF button.
2. Press the TEST button until the arrow indicates the desired test-method memory.
3. Compare the first three digits of the batch number of the Reflectoquant® pack. This reference number must be the same as the one displayed (otherwise recalibrate, see page 10).
4. Press the START button. The waiting time appears in the display, e.g. 120 sec (A). Proceed as described in the respective Reflectoquant® test package insert. Press the START button again to activate the stopwatch function (B). When the waiting time has expired, a strip symbol appears in the display and the reaction time of the test, e.g. 60 sec, is shown (C).
5. Immerse the test strip into your sample as described in the Reflectoquant® package insert (Reflectoquant® pack) and at the same time press the START button (A). Carefully allow excess liquid to run off via the long edge of the strip onto an absorbent paper towel. The countdown function of the instrument is activated, the remaining reaction time is shown (countdown) (B).
6. An acoustic signal and the blinking strip five seconds before the end of the reaction time requests you to insert the strip all the way into the strip adaptor. The strip can also be inserted before the five seconds have expired (see the corresponding package insert).
7. The measurement result is shown in the display in the units for the respective test, e.g. mg/l, and is automatically stored.






Measurement procedure C

This procedure is used to evaluate cells. Insert the cell adapter.

How?

1. Switch on the instrument by pressing the ON/OFF button.
2. Press the TEST button until the arrow indicates the desired test-method memory.
3. Compare the first three digits of the batch number of the Reflectoquant® pack. This reference number must be the same as the one displayed (otherwise recalibrate, see page 10).
4. Press the START button. The lid must be closed. Where applicable, the waiting time appears in the display, e.g. 240 seconds.
5. Proceed as described in the respective Reflectoquant® plus test package insert. Press the START button anew to activate the stopwatch function (A). When the waiting time has expired, the reaction time appears in the display, e.g. 5 seconds (B).
6. Insert the cell containing the blank solution into the cell adapter, close the lid, and press the START button (A). An acoustic signal sounds while the blank is being measured.
Two alternately blinking cell symbols on the display then prompt you to insert the cell containing the sample solution (B).
7. After inserting the cell, close the lid and press the START button anew (A). An acoustic signal sounds while the sample is being measured.
The measurement result is shown in the display in the units for the respective test, e.g. mg/l, and is automatically stored (B).
8. After measurement, immediately remove the cell from the cell adapter.

Important

Interrupting the measurement procedure (menu guidance), to perform measurements directly without the countdown function of the instrument.

When?

1. To repeat a measurement e.g. when the test strip / cell was inserted incorrectly or too late.
2. For serial measurements.

How?

Once the first measurement has been completed, a further measurement can be carried out simply by pressing the START button. In all these cases the instrument makes the measurements immediately.

The result is displayed and stored automatically.

Remarks regarding serial measurements

If you wish to measure e.g. several nitrate samples, it is advisable to follow the procedure below. For repeat serial measurement the countdown function is not available and an additional stopwatch is needed.

1. Run the standard measurement (procedure A) once.
2. Immerse separate test strips into your sample at 15-second intervals. Carefully allow excess liquid to run off via the long edge of the strip onto an absorbent paper towel and allow each strip to react outside the instrument.
3. After completion of the reaction time (60 sec) insert each strip into the instrument. Press the START button to perform the measurement (i.e. a strip is inserted and measured every 15 seconds).

Attention

The reaction times of the tests must be monitored precisely. Differences from the default value may cause (small) errors in the result due to changes in strip color with time.

The instrument switches off after 2 minutes if no key is pressed within this time period.

5 Handling measurement results

When?

1. To store measurement results.
2. To recall measurement results.
3. To transfer measurement results to a PC.
4. To delete stored measurement results.

Storing of measurement results

How?

Measurement results are automatically stored together with the time and date.

The instrument stores up to 50 results.

When the batch number is shown the number of available memory slots, e.g. F:21, is displayed instead of the time and date.

Important

The display of F:00 indicates that all 50 storage slots are occupied. The next measurement result overwrites the oldest one.


Displaying stored data

How?

1. Switch on the instrument by pressing the ON/OFF button.
2. Press the TEST button until the arrow indicates the desired test-method memory. The reference code is shown.

3. Press the MEM button. The last stored result of the method in question is shown. If no data are stored, an acoustic signal sounds. The display remains unchanged.
4. Repeated pressing of the MEM button displays the stored data in reverse time and date order (A). The display of the reference code in conjunction with the sound of the acoustic signal indicates that no further data are stored (B).
5. To select stored data of another test-method memory press the TEST button to set the pointer. Then press the MEM key to display the data.
6. To finish the operation, press the ON/OFF (or the TEST) button.


2


4B


3

5

4A

Transferring data to a PC

When?

For quality assurance and/or documentation of measurement results.

How?

The instrument is connected to the PC via a dedicated interface.

You also require the RQdata package, consisting of software and a data-transfer cable, Cat. No. 1.16998.0001.

Deleting data

When?

1. To delete the last result.
2. In the case of an incorrect measurement.
3. To delete all measurement results.

How?

When the reference code is displayed, press the MEM button for three seconds. This deletes the most recent stored result. During this procedure, the result to be deleted blinks in the display and an acoustic signal sounds.

Deletion is complete when the display shows 000. Press the TEST button.

The reference code appears again.

Important


When recalibrating with the bar code, all measurement results stored in the test-method memory are deleted!

6 Maintenance of the instrument

To obtain consistently exact measurement results please follow these instructions for cleaning and maintenance.

Handling

1. Please treat this instrument with the same care as you do with all other electronic devices. Take care to ensure that no liquid enters the instrument case, since this may result in damage that is not covered by the warranty.
2. Do not expose the instrument to excessively humid conditions or to excessive heat or cold (see the technical data).
3. If the external case of the measurement instrument needs cleaning, please wipe the case and the display carefully with a moist cloth.
4. Clean the optics, the strip adaptor, and the cell adaptor only with distilled water or ethanol.

Cleaning the strip adapter


When?

The strip adapter must be cleaned immediately if an error message (OPT or ERR) is displayed.

The strip adapter should be thoroughly cleaned (dismantling the strip adapter into its individual components) at the end of each working day.

How?

1. Switch off the instrument by pressing the ON/OFF button.
2. Carefully pull the strip adapter out of the case.
3. Dismantle the adapter into its three components.
4. Clean the components with water and a mild detergent. If necessary use ethanol. Never clean the internal standard (light grey plastic component) with abrasive cleaner!
5. Dry the components carefully and reassemble the adapter.
6. Reinsert the strip adapter into the instrument.

Cleaning the cell adapter

When?

The strip adapter must be cleaned immediately if any sample solution is inadvertently spilled in the measurement chamber.

How?

1. Switch off the instrument by pressing the ON/OFF button.
2. Carefully pull the cell adapter out of the case.
3. Clean the components with water and a mild detergent.
5. Dry the cell adapter and optics immediately.
6. Reinsert the cell adapter into the instrument.

Recalibration

The recalibration set consists of an internal standard (light grey plastic component), a bar-code strip for calibration, and a white calibration strip.


When?

1. Once a week as regular maintenance.
2. In the case of suspect measurement results.
3. After a change of the strip adapter and/or the internal standard or batteries.
4. After severe mechanical distress (e.g. fall).
5. In the case of the error message (OPT or ERR).

How?

1. Clean the strip adapter thoroughly. Take care that the internal standard has not changed color (if necessary exchange the internal standard, Cat. No. 1.16953.0001).
2. Keep the instrument at room temperature for at least 30 minutes.
3. Insert the adapter and switch on the instrument.
4. Insert the bar-code strip for calibration. CAL appears in the display.

5. Insert the calibration strip of the recalibration set into the strip adapter as shown in the diagram, and ensure that it is correctly positioned. Take utmost care to ensure that the calibration strip is clean.
6. Press the START key. The instrument is recalibrated, the CAL message disappears from the display, and the instrument switches off automatically.
7. Store the recalibration set in the supplied packaging materials, taking care to protect it from light.

Important

When the batteries are low (display) the instrument will not accept the bar-code strip for calibration. Change the batteries first. If the calibration strip is not inserted into the measurement chamber, an acoustic signal sounds after the START button is pressed. The blinking CAL symbol in the display prompts you to insert the calibration strip.

7 Trouble-shooting

The following table shows explanations of the error messages and tips on how to avoid incorrect measurements.

The most commonly encountered problems are caused by:

- The test strips or cells not being inserted correctly.
- Cells being dirty or strongly scratched.
- The reaction time of the strip not being right.
- Incorrect use of the tests (e.g. strips / tests not properly stored, pH range not correctly adjusted).

Problem	Cause	Solution
• No display	<ul style="list-style-type: none"> • Batteries not inserted properly • Batteries are drained 	<ul style="list-style-type: none"> • Insert the batteries correctly • Use new batteries
• Display LO	• Concentration is lower than the measurement range of the test	• State the result is smaller than the lowest measurement value
• Display HI	• Concentration is higher than the measurement range of the test	• Dilute the sample; do not forget to consider the factor for dilution when reporting the result!
• Display - - -	• No calibration programmed	• Calibrate instrument

Problem	Cause	Solution
• Display OPT	• Optics or strip adapter is dirty	• Clean optics and adapter and switch the RQflex® plus 10 on and then off again; recalibrate, see page 26
• Display ERR	<ul style="list-style-type: none"> • Optical error • Strip adapter is assembled incorrectly 	• Clean optics and adapter; recalibrate instrument, see page 26; if necessary, call customer service
• Display E-1	<ul style="list-style-type: none"> • Surrounding is too bright • Strip adapter is assembled incorrectly 	• Repeat measurement in a darker place; if necessary, close the lid of the instrument before switching on
• Display E-2	• Result is not plausible	• Repeat measurement

8 Service/Warranty

Our instruments are 100% quality-controlled, i.e. each and every individual instrument has been tested before leaving the factory. The RQflex® plus 10 is maintenance-free. The detailed instructions should allow for problem-free operation.

In the event that you are nevertheless faced with an unsolvable problem, please contact the technical customer service of our supplier. There you will receive all the information necessary to help you solve the problem.

Important

In the event that you wish to return the instrument, whatever the reason, please contact our technical customer service first.

Warranty

1. The manufacturer grants for this RQflex® plus 10 a total replacement warranty of 24 months from the date of purchase.
2. All warranty work is provided by our customer service unit. In the event of a proven production or material defect, we shall provide you with a new or as-new replacement RQflex® plus 10 free of charge.
3. Improper handling results in the warranty becoming invalid. No water or any other liquid must be allowed to enter the instrument. In such cases the warranty becomes invalid.
4. All warranty rights become invalid in the event that the purchaser or an unauthorized third person repairs or opens the instrument or if changes are made to the warranty slip. Only the customer-service unit is authorized to service components.
5. If the instrument is sent in for repair under the terms of the warranty, a copy of the invoice or of another proof of purchase must be enclosed in all cases.

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EMD Millipore Corporation, 290 Concord Road, Billerica,